

COVID-19 Protection Equipment Product Catalogue

**Protect your self and your loved ones by using our wide
range of sanitization products**

**Our range of products include sanitization devices for
home and office use and personal protective gear**

Available at best prices

Retail and Bulk wholesale purchase options

For enquiries

call us on 8850489889 or email at sales@arshderma.com

Handheld Contactless Thermal Scanner

360 Degree View



Call for best price

Speical bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com

Fog Machine for Home and Small Office



Call for best price

Speical bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com

Personal Protective Equipment



PPE Kit
 PERSONAL PROTECTION EQUIPMENT
 FOR CORONA VIRUS

NO	DESCRIPTION	QTY
1	Cover All Suit with Hoodie (90 GSM)	1 No.
2	3 ply Mask	1 No.
3	Hand Gloves	1 Pair
4	Face Shield	1 No.
5	Shoe Cover	1 Pair
6	Bio Waste Bag	1 No.

Call for best price

Special bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com

Non – Medical N95 Mask, Face Shield and Normal Mask

KN95
 Non-Medical
 Ear Hanging Stero Mask Professional
 Protection And Health Care

CE

FUNCTION
 This product is designed to provide protection against the spread of respiratory droplets.

USAGE
 1. Before using the mask, clean your hands with soap and water for at least 20 seconds.
 2. Hold the mask by the top edge and place it over your nose and mouth.
 3. Press the bridge of the mask to fit snugly against your face.
 4. Do not touch the front of the mask.
 5. After use, discard the mask in a bin or wash it with soap and water for at least 20 seconds.

PRODUCT INFORMATION
 PRODUCT NAME: PROTECTIVE MASK
 PRODUCT FEATURE AND COMPOSITION: THIS PRODUCT IS LABOURATED FROM THE BEST QUALITY OF THE MARKET.
 SIZE: ADULT
 REFILL/RECHARGE/RECHARGE OR REUSE INFO:
 REUSABLE PRODUCT. SEE USER BOX CERTIFICATE.
 PACKING & STORAGE:
 PACKING: STERILE IN A VACUUM, BPA FREE AND NON-CORROSIVE GAS PLACE.
 STORAGE: STORE IN A COOL, DRY PLACE. DO NOT EXPOSE TO DIRECT SUNLIGHT.
 PLEASE USE WITH CAUTION.
 THIS PRODUCT CANNOT BE USED AS A GAS MASK IN A TOXIC ENVIRONMENT.

20PCS

1. Non-Woven(50g)
 2. Melt-Blown (95/25g)
 3. Hot-Air Cotton
 4. Melt-Blown (95/25g)
 5. Non-Woven(25g)

Magiclayers **Face Shield** **CE**

How to use ?
 1. Take out the plastic head band and transparent sheet from the box.
 2. Align four holes of the sheet with the four notches on the head band and press so that the sheet snaps into the notches.

Features :
 - CE Certified
 - Full face coverage (Eyes, Nose, Mouth)
 - No need of elastic or string, fits directly on the head
 - Light weight and comfortable
 - Reusable

Cleaning methods : Use any one method
 Separate the headband and sheet to undertake cleaning process.

Wash with soap water for 20 secs
Clean with 70% alcohol
(Disinfect Only)
 Autoclave the headband and replace the sheet

250 microns PETG Sheet

ABS Headband

Call for best price

Special bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com

VENUS – N95 Masks

New! VENUS RPD714N95



Call for best price

Special bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com

IMMUNITY BOOSTER & SKIN CARE



Call for best price

Special bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com

Contact-less Sanitization Tunnel for Offices



Call for best price

Speical bulk discount available based on purchase quantity

For enquiries

call us on 8850489889 or email at sales@arshderma.com